

	<a href="http://www.e-labeling.eu/CYL10508">www.e-labeling.eu/CYL10508</a>														
	+800 135 79 135														
8:00 - 17:00 GMT+1	GR 00 800 161 220 577 99														
<table border="1"> <tr> <td>M</td><td>T</td><td>W</td><td>T</td><td>F</td><td>S</td><td>S</td> </tr> <tr> <td>☒</td><td>☒</td><td>☒</td><td>☒</td><td>☒</td><td>☒</td><td>☒</td> </tr> </table>		M	T	W	T	F	S	S	☒	☒	☒	☒	☒	☒	☒
M	T	W	T	F	S	S									
☒	☒	☒	☒	☒	☒	☒									

## Magnet Tray

The magnet tray is used in combination with Immuknow® – the Cylex™ Immune Cell Function Assay. For additional information on the correct use of the magnet tray, please see the Immuknow package insert.

### The Magnet Tray Consists of 2 Components:

1. The strip holder (top portion) accommodates up to six 1 x 8 strips from the assay plate supplied in the assay kit. This portion is completely separable from the bottom in order to remove samples from the magnetic field. There is only one recommended way of removing the 1 x 8 strips from the assay plate and placing them into the strip holder so that both ends of the strip are properly and completely seated.
2. The magnet base (lower portion) contains 3 rows of strong magnets. It is designed to attract magnetic particles to the sides of each well of the assay plate within a short time to facilitate washing the wells.

### Cleaning Precautions:

1. Do not disinfect the magnet tray using alcohol or an autoclave.
2. Cleaning with other disinfectants (e.g., 10% bleach) should be followed by thorough rinsing with water to remove residue that may cause corrosion.
3. Dry the magnet tray completely before storage.

**Caution: Keep the magnet tray away from magnet-sensitive materials (e.g., magnetic recording material, computer disks, credit cards).**



Qarad b.v.b.a.  
Volmolenheide 13  
B-2400 Mol, Belgium



Be in the Know.

**Immuknow®**